



SAMM50 Antibody

Product Code	CSB-PA020692GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y512
Immunogen	Human SAMM50
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	sorting and assembly machinery component 50 homolog (S. cerevisiae);SAMM50;CGI-51;FLJ35825;FLJ42905;FLJ99036;OMP85;SAM50;T OB55;TRG-3;YNL026W ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SAMM50
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.