



PPP6R3 Antibody

Product Code	CSB-PA020701GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5H9R7
Immunogen	Human SAPS3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SAPS domain family, member 3;SAPS3;C11orf23;DKFZp781E17107;DKFZp781E2374;DKFZp781O2362;FLJ11058;FLJ43065;KIAA1558;MGC125711;MGC125712;PP6R3;SAP190;SAPL;SAPLa ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP6R3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.