



SCFD1 Antibody

Product Code	CSB-PA020806GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVM8
Immunogen	Human SCFD1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	sec1 family domain containing 1;SCFD1;C14orf163;KIAA0917;RA410;SLY1;STXBP1L2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SCFD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.