



SCG5 Antibody

Product Code	CSB-PA020810GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05408
Immunogen	Human SCG5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	secretogranin V (7B2 protein);SCG5;7B2;P7B2;SGNE1;SgV ;SECRETOGRANIN V
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SCG5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.