



SERBP1 Antibody

Product Code	CSB-PA021038GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NC51
Immunogen	Human SERBP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	SERPINE1 mRNA binding protein 1;SERBP1;CGI-55;CHD3IP;DKFZp564M2423;FLJ90489;HABP4L;PAI- RBP1;PAIRBP1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SERBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.