



SERINC1 Antibody

Product Code	CSB-PA021046GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRX5
Immunogen	Human SERINC1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	serine incorporator 1; SERINC1; KIAA1253; TDE1L; TDE2; TMS-2 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SERINC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.