



SERP1 Antibody

Product Code	CSB-PA021051GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6X1
Immunogen	Human SERP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	stress-associated endoplasmic reticulum protein 1;SERP1;MGC117327;MGC133321;MGC133322;RAMP4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SERP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.