



SERPINA5 Antibody

Product Code	CSB-PA021061GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05154
Immunogen	Human SERPINA5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 5;SERPINA5;PAI3;PCI;PLANH3;PROCI ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SERPINA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.