



# SERPINH1 Antibody

<b>Product Code</b>	CSB-PA021087GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P50454
<b>Immunogen</b>	Human SERPINH1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding protein 1);SERPINH1;AsTP3;CBP1;CBP2;HSP47;PIG14;PPROM;RA-A47;SERPINH2;gp46 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	SERPINH1
<b>Target Details</b>	This gene encodes a member of the serpin superfamily of serine proteinase inhibitors. Its expression is induced by heat shock. The protein localizes to the endoplasmic reticulum lumen and binds collagen; thus it is thought to be a molecular chaperone involved in the maturation of collagen molecules. Autoantibodies to this protein have been found in patients with rheumatoid arthritis.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.