



SETD1B Antibody

Product Code	CSB-PA021102GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPS6
Immunogen	Human SETD1B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SET domain containing 1B;SETD1B;FLJ20803;KIAA1076;KMT2G;Set1B ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SETD1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.