



SFTPb Antibody

Product Code	CSB-PA021173GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07988
Immunogen	Human SFTPb
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Form	Liquid
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	surfactant protein B; SFTPb; PSP-B; SFTB3; SFTP3; SMDP1; SP-B ;
Immunogen Species	Homo sapiens (Human)
Target Names	SFTPb
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.