



SHANK1 Antibody

Product Code	CSB-PA021247GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y566
Immunogen	Human SHANK1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	0.1M NaHCO ₃ ,0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SH3 and multiple ankyrin repeat domains 1;SHANK1;SPANK-1;SSTRIP;synamon ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SHANK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.