



SIGLEC5 Antibody

Product Code	CSB-PA021299GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15389
Immunogen	Human SIGLEC5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	sialic acid binding Ig-like lectin 5; SIGLEC5; CD170; CD33L2; OB-BP2; OBBP2; SIGLEC-5 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SIGLEC5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.