



SIRT1 Antibody

Product Code	CSB-PA021339GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EB6
Immunogen	Human SIRT1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF, IP; Recommended dilution: WB:1:500-1:1000, IHC:1:50-1:200, IF:1:50-1:500
Form	Liquid
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	sirtuin (silent mating type information regulation 2 homolog) 1 (S. cerevisiae);SIRT1;SIR2L1 ;
Immunogen Species	Homo sapiens (Human)
Target Names	SIRT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.