



SIX2 Antibody

Product Code	CSB-PA021349GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPC8
Immunogen	Human SIX2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	SIX homeobox 2;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SIX2
Target Details	This gene is a member of the vertebrate gene family which encode proteins homologous to the Drosophila sine oculis homeobox protein. The encoded protein is a transcription factor which, like other members of this gene family, may be involved in limb or eye development.