



SKAP2 Antibody

Product Code	CSB-PA021355GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75563
Immunogen	Human SKAP2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	src kinase associated phosphoprotein 2;SKAP2;MGC10411;MGC33304;PRAP;RA70;SAPS;SCAP2;SKAP-HOM;SKAP55R ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SKAP2
Target Details	This protein belongs to the src family kinases. This protein is similar to the src kinase associated phosphoprotein 1. It is an adaptor protein that is thought to play an essential role in the src signaling pathway in various cells. It inhibits PTK2B/RAFTK activity and regulates alpha-synuclein phosphorylation.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.