



SLC20A1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA021439GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8WUM9 |
| Immunogen | Human SLC20A1 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | solute carrier family 20 (phosphate transporter), member 1;SLC20A1;DKFZp686J2397;FLJ41426;GLVR1;Glvr-1;PIT1;PiT-1 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | SLC20A1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |