



SLC8B1 Antibody

Product Code	CSB-PA021474GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6J4K2
Immunogen	Human SLC24A6
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 24 (sodium/potassium/calcium exchanger), member 6; SLC24A6; FLJ22233; NCKX6; NCLX ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC8B1