



SLC25A25 Antibody

Product Code	CSB-PA021494GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6KCM7
Immunogen	Human SLC25A25
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 25 (mitochondrial carrier; phosphate carrier), member 25;SLC25A25;KIAA1896;MCSC;MGC105138;MGC119514;MGC119515;MGC119516;MGC119517;PCSC;RP11-395P17.4;SCAMC-2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC25A25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.