



SLC27A3 Antibody

Product Code	CSB-PA021535GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5K4L6
Immunogen	Human SLC27A3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	solute carrier family 27 (fatty acid transporter), member 3;SLC27A3;ACSVL3;FATP3;MGC4365;VLCS-3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC27A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.