



SLC38A4 Antibody

Product Code	CSB-PA021619GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969I6
Immunogen	Human SLC38A4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 38, member 4;SLC38A4;ATA3;FLJ10191;MGC126876;NAT3;PAAT ;SNAT4;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC38A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.