



SLC44A1 Antibody

Product Code	CSB-PA021648GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWI5
Immunogen	Human SLC44A1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 44, member 1;SLC44A1;CD92;CDW92;CHTL1;CTL1;RP11-287A8.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC44A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.