



SLC7A1 Antibody

Product Code	CSB-PA021709GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30825
Immunogen	Human SLC7A1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 7 (cationic amino acid transporter, y+ system), member 1;SLC7A1;ATRC1;CAT-1;ERR;HCAT1;REC1L ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC7A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.