



SLC9A8 Antibody

Product Code	CSB-PA021739GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2E8
Immunogen	Human SLC9A8
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 9 (sodium/hydrogen exchanger), member 8;SLC9A8;DKFZp686C03237;FLJ42500;KIAA0939;MGC138418;NHE8 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC9A8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.