



SLK Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA021776GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9H2G2 |
| Immunogen | Human SLK |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA,WB,IF |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | STE20-like kinase (yeast);SLK;KIAA0204;LOSK;MGC133067;STK2;bA16H23.1;se20-9 ; |
| Product Type | Purified Rabbit Anti Human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | SLK |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |