



SMAD5 Antibody

Product Code	CSB-PA021790GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99717
Immunogen	Human SMAD5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	SMAD family member 5;SMAD5;DKFZp781C1895;DKFZp781O1323;Dwfc;JV5-1;MADH5 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SMAD5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.