



MIEF1 Antibody

Product Code	CSB-PA021827GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQG6
Immunogen	Human SMCR7L
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Smith-Magenis syndrome chromosome region, candidate 7-like;SMCR7L;FLJ20232;HSU79252;dJ1104E15.3 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MIEF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.