



SMS Antibody

Product Code	CSB-PA021854GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52788
Immunogen	Human SMS
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	spermine synthase;SMS;MRSR;SPMSY;SRS;SpS ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SMS
Target Details	This protein belongs to the spermidine/spermine synthases family. This gene encodes an ubiquitous enzyme of polyamine metabolism.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.