



SMU1 Antibody

Product Code	CSB-PA021858GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2TAY7
Immunogen	Human SMU1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	smu-1 suppressor of mec-8 and unc-52 homolog (C. elegans);SMU1;BWD;DKFZp761L0916;FLJ10805;FLJ10870;FLJ11970;MGC117363;RP11-54K16.3;SMU-1;fSAP57 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SMU1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.