



SMYD2 Antibody

Product Code	CSB-PA021863GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRG4
Immunogen	Human SMYD2
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SET and MYND domain containing 2; SMYD2; HSKM-B; KMT3C; MGC119305; ZMYND14 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SMYD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.