



SP100 Antibody

Product Code	CSB-PA022439GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23497
Immunogen	Human SP100
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	SP100 nuclear antigen; SP100; DKFZp686E07254; FLJ00340; FLJ34579 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SP100
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.