



SPAG16 Antibody

Product Code	CSB-PA022462GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0X2
Immunogen	Human SPAG16
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	sperm associated antigen 16;SPAG16;DKFZp666P1710;FLJ22724;FLJ37717;MGC87036;PF20;WDR29 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SPAG16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.