



SRC Antibody

Product Code	CSB-PA022650GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12931
Immunogen	Human SRC
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IP,IF;WB:1:500-1:2000,IP:0.5-4.0ug,IHC:1:100-1:400,IF:1:10-1:100
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	v-src sarcoma (Schmidt-Ruppin A-2) viral oncogene homolog (avian);SRC;ASV;SRC1;c-SRC;p60-Src ;
Immunogen Species	Homo sapiens (Human)
Target Names	SRC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.