



SRPR Antibody

Product Code	CSB-PA022686GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08240
Immunogen	Human SRPR
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	signal recognition particle receptor (docking protein);SRPR;DP;MGC17355;MGC3650;MGC9571;SRP-alpha;Sralpha ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SRPRA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.