



SSSCA1 Antibody

Product Code	CSB-PA022722GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60232
Immunogen	Human SSSCA1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	Sjogren syndrome/scleroderma autoantigen 1;SSSCA1;p27 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SSSCA1
Target Details	This antigen is recognized by a subset of anti-centromere antibodies from patients with scleroderma and/or Sjogren s syndrome. Subcellular localization has not yet been established.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.