



ST8SIA3 Antibody

Product Code	CSB-PA022776GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43173
Immunogen	Human ST8SIA3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ST8 alpha-N-acetyl-neuraminide alpha-2,8-sialyltransferase 3;ST8SIA3;SIAT8C;ST8SiaIII ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ST8SIA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.