



STBD1 Antibody

Product Code	CSB-PA022820GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95210
Immunogen	Human STBD1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	starch binding domain 1; STBD1; FLJ41801; GENEX3414; GENX-3414; Genethonin 1
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STBD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.