



STEAP4 Antibody

Product Code	CSB-PA022826GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q687X5
Immunogen	Human STEAP4
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	STEAP family member 4; STEAP4; DKFZp666D049; FLJ23153; STAMP2; TIARP; TNFAIP9 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STEAP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.