



STK11IP Antibody

Product Code	CSB-PA022834GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N1F8
Immunogen	Human STK11IP
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	serine/threonine kinase 11 interacting protein;STK11IP;KIAA1898;LIP1;LKB1IP;STK11IP1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STK11IP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.