



STXBP2 Antibody

Product Code	CSB-PA022906GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15833
Immunogen	Human STXBP2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	syntaxin binding protein 2;STXBP2;Hunc18b;MUNC18-2;UNC18-2;UNC18B;pp10122 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STXBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.