



SULT1B1 Antibody

Product Code	CSB-PA022937GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43704
Immunogen	Human SULT1B1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	sulfotransferase family, cytosolic, 1B, member 1; SULT1B1; MGC13356; ST1B2; SULT1B2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SULT1B1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.