



SURF2 Antibody

Product Code	CSB-PA022964GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15527
Immunogen	Human SURF2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	surfeit 2;SURF2;SURF-2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SURF2
Target Details	This gene is located in the surfeit gene cluster, a group of very tightly linked genes that do not share sequence similarity. The gene shares a bidirectional promoter with SURF1, which is located on the opposite strand. The function of the encoded protein from this housekeeping gene has not yet been determined.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.