



# SURF2 Antibody

<b>Product Code</b>	CSB-PA022964GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q15527
<b>Immunogen</b>	Human SURF2
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	surfeit 2;SURF2;SURF-2 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	SURF2
<b>Target Details</b>	This gene is located in the surfeit gene cluster, a group of very tightly linked genes that do not share sequence similarity. The gene shares a bidirectional promoter with SURF1, which is located on the opposite strand. The function of the encoded protein from this housekeeping gene has not yet been determined.