



SV2C Antibody

Product Code	CSB-PA022980GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q496J9
Immunogen	Human SV2C
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	synaptic vesicle glycoprotein 2C;SV2C;MGC118874;MGC118875;MGC118876;MGC118877 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SV2C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.