



# TBC1D22A Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA023196GA01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                        |
| <b>Uniprot No.</b>         | Q8WUA7   |
| <b>Immunogen</b>           | Human TBC1D22A   |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human,Mouse,Rat  |
| <b>Tested Applications</b> | ELISA,WB,IHC   |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| <b>Purification Method</b> | Antigen Affinity purified  |
| <b>Isotype</b>             | IgG  |
| <b>Alias</b>               | TBC1 domain family, member 22A;TBC1D22A;C22orf4;HSC79E021 ;                          |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody                                       |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | TBC1D22A   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.          |