



TCF7L1 Antibody

Product Code	CSB-PA023306GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCS4
Immunogen	Human TCF7L1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	transcription factor 7-like 1 (T-cell specific, HMG-box);TCF7L1;TCF-3;TCF3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TCF7L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.