



TCTN1 Antibody

Product Code	CSB-PA023338GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2MV58
Immunogen	Human TCTN1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	tectonic family member 1;TCTN1;FLJ21127;TECT1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TCTN1
Target Details	This gene encodes a member of the tectonic family of secreted and transmembrane proteins. The orthologous gene in mouse is required for formation of most ventral cell types. It functions downstream of smoothed and rab23 to modulate hedgehog signal transduction. Multiple transcript variants encoding different isoforms have been found for this gene.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.