



TEX12 Antibody

Product Code	CSB-PA023402GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXU0
Immunogen	Human TEX12
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	testis expressed 12;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TEX12
Target Details	This gene is similar to a mouse gene that is expressed in the testis.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.