



TFDP2 Antibody

Product Code	CSB-PA023426GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14188
Immunogen	Human TFDP2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	transcription factor Dp-2 (E2F dimerization partner 2);TFDP2;DP2;Dp-2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TFDP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.