



TFEB Antibody

Product Code	CSB-PA023429GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19484
Immunogen	Human TFEB
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IP,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	transcription factor EB;TFEB;AlphaTFEB;TCFEB;bHLHe35 ;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TFEB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.